

## Protocol for Detection of Antibody Internalization

### 1. Incubate pH Sensitive IgG Labeling Reagent with Antibody

- 1.1 Dissolve pH Sensitive IgG labeling reagent (22310) in freeze-dried powder according to COA using deionized water.
- 1.2 Mix antibody with the pH Sensitive IgG labeling reagent (22310) in 9:1 ratio by weigh (2:1 ratio by molarity), and then dilute the mixture of the antibody and pH Sensitive IgG labeling reagent to two-fold of a designed concentration using complete culture medium.
- 1.3 Incubate at room temperature in dark for 1 hour to obtain antibody-pH Sensitive IgG labeling reagent (22310) complex.

### 2. Incubate the Complex of Antibody-pH Sensitive IgG Labeling Reagent with Cells

- 2.1 Collect, wash cells and adjust cell concentration using complete culture medium ( $1-2 \times 10^5$  cells/mL for suspend cells or  $0.5-1 \times 10^5$  cells/mL for adherent cells). Add 100  $\mu$ L cell suspension to each hole of the 96 well plate.
- 2.2 Add 100  $\mu$ L Antibody-pH Sensitive IgG labeling reagent (22310) complex to each hole. Then incubate in a 5% CO<sub>2</sub> incubator at 37 °C.

### 3. Detect the Internalization with Flow Cytometry

After incubation for 18-24 hours, the internalization effect of the antibody was detected by flow cytometry using FITC or AF488 channels.

---

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS