

## **Product Description**

Pioneering GTPase and Oncogene Product Development since 2010

## SARS-COV-2 (2019-NCOV) NUCLEOCAPSID, N-HIS TAG

货号: 11271

产品全名: SARS-CoV-2 (2019-NCoV) Nucleocapsid

规格: 10/50/100 μg

基因符号 Nucleocapsid protein;P;rotein N

目标蛋白: Nucleocapsid UNIPROT ID: PODTC9

描述: "SARS-CoV-2 (2019-nCoV) Nucleocapsid

背景: Coronavirus contain most of nucleocapsid protein. Coronavirus nucleoproteins (N proteins) localize to the cytoplasm and the nucleolus, a subnuclear structure, in both virus-infected primary cells and in cells transfected with plasmids that express N protein. The nucleolus is the site of ribosome biogenesis and sequesters cell cycle regulatory complexes. Two of the major components of the nucleolus are fibrillarin and nucleolin. These proteins are involved in nucleolar assembly and ribosome biogenesis and act as chaperones for the import of proteins into the nucleolus. Regarding of the conservation of N protein sequence and its strong immunogenicity, the N protein of coronavirus is a tool for diagnostic.

物种/宿主: E.coli

分子量: The protein has a predicted molecular mass of 49.4 kDa after removal of the signal peptide.

分子特征: 6\*His Tag Nucleocapsid protein(Met1-Ala419)

纯化: The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCI, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as

protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.





Figure 1. SARS-CoV-2 (2019-nCoV) Nucleocapsid Protein, His Tag on SDS-PAGE under reducing condition.

## SARS-CoV-2 (2019-nCoV) Nucleocapsid, His Tagged protein ELISA

 $0.2~\mu g$  of SARS-CoV-2 Nucleocapsid protein, His Tagged per well



Concentration of Anti-SARS-CoV-2 Nucleocapsid Antibody (ng/ml)

Figure 2. ELISA plate pre-coated by 2  $\mu g/ml$  (100  $\mu l/well$ ) SARS-CoV-2 (2019nCov) Nucleocapsid, His Tag protein (1127) can bind Anti-SARS-CoV-2 Nucleocapsid Antibody 28304], [getskuurl sku=28305 in a linear range of 0.122-15.625 ng/ml.