

## **Product Description**

Pioneering GTPase and Oncogene Product Development since 2010

## SARS-COV-2 (2019-NCOV) S1 PROTEIN NTD, HFC TAG

货号: 11299

产品全名: SARS-CoV-2 (2019-NCoV) SI Protein NTD

规格: 10/50/100 µg

基因符号 SI protein NTD;pike protein SI NTD;etaCoV SI-NTD;OVID-19

目标蛋白: S1 protein NTD UNIPROT ID: P0DTC2

描述: Recombinant SARS-CoV-2 (2019-nCoV) SI protein NTD with C-terminal human Fc tag

背景: SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 59.0 kDa after removal of the signal peptide. The apparent molecular mass of S1-NTD-hFc is approximately 70-100 kDa due to glycosylation.

分子特征: S1 protein NTD(Ser13-Leu303) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at −20°C to −80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at −80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.





Figure 1. SARS-CoV-2 (2019-nCoV) S1 protein NTD, hFc Tag on SDS-PAGE under reducing condition.



SARS-CoV-2 (2019-nCoV) S1 NTD, hFc Tagged protein(µg/ml)

Figure 2. ELISA plate pre-coated by 2 μg/ml (100 μl/well) Human AXL, His tagged protein 11193 can bind SARS-CoV-2 (2019-nCoV) S1 protein NTD, hFc Tagged protein "ΜΔ11299" in a linear range of 0.0012-0.16 μg/ml.