

## SARS-COV-2 (BETA) S PROTEIN RBD, HFC TAG

货号: 11339

产品全名: SARS-CoV-2 (Beta) S Protein RBD

规格: 10/50/100 µg

基因符号 SARS-CoV-2 B.1.351 (Beta) Spike RBD 蛋白

目标蛋白: S protein RBD

**UNIPROT ID:** P0DTC2

**描述:** "Recombinant SARS-CoV-2 S protein RBD(K417N

**背景:** SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 51.3 kDa after removal of the signal peptide. The apparent molecular mass of RBD(K417N, E484K, N501Y)-hFc is approximately 55-70 kDa due to glycosylation.

**分子特征:** S protein RBD(K417N, E484K, N501Y)(Arg319-Phe541) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 90% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. SARS-CoV-2 (2019-nCoV) S protein RBD (K417N, E484K, N501Y), hFc Tag on SDS-PAGE under reducing condition.

**CoV-2 (Beta) S-RBD, hFc Tagged protein ELISA**  
0.2 µg of CoV-2 (Beta) S-RBD, hFc tagged protein per well




Figure 2. ELISA plate pre-coated by 2 µg/ml (100 µl/well) SARS-CoV-2 (Beta) S protein RBD, hFc Tag (11339) can bind Human ACE2 Protein, mFc Tag (11195) in a linear range of 12.346–3000 ng/mL.

**CoV-2 (Beta) S-RBD, hFc Tagged protein ELISA**  
0.2 µg of CoV-2 (Beta) S-RBD, hFc tagged protein per well




Figure 3. ELISA plate pre-coated by 2 µg/ml (100 µl/well) SARS-CoV-2 (Beta) S protein RBD, hFc Tag (11339) can bind Anti-SARS-CoV-2 RBD antibody (DM55), Rabbit mAb (28354) in a linear range of 0.244–15.625 ng/mL.