

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

SARS-COV-2 (OMICRON BA.5) S PROTEIN RBD, HFC TAG

货号: 11810

产品全名: SARS-CoV-2 (Omicron BA.5) S Protein RBD

规格: 10/50/100 µg

基因符号 SARS-CoV-2 BA.5 (Omicron) Spike RBD 蛋白

目标蛋白: S protein RBD UNIPROT ID: PODTC2

描述: Recombinant SARS-CoV-2 (Omicron BA.5) S-RBD protein with C-terminal human Fc tag

背景: SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 51.4 kDa after removal of the signal peptide. The apparent molecular mass of S-RBD(Omicron BA.5)-hFc is approximately 55-70 kDa due to glycosylation. 分子特征: S protein RBD(G339D, S371F, S373P, S375F, T376A, D405N, R408S, K417N,N440K, L452R,S477N, T478K, E484A,F486V, Q498R, N501Y,

Y505H)(Arg319-Phe541) hFc(Glu99-Alg330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010




Figure 1. SARS-CoV-2 (Omicron BA.5) S protein RBD, hFc Tag on SDS-PAGE under reducing condition.